

EasyStand Zing Prone Support Package Feature Descriptions

FEATURE	IMAGE	DESCRIPTION	MOD	MAX
Multi-Adjustable Footplates with Foot Holders		Foot plates are adjustable forward/aft +/75" (2cm) +extra 1.5" (3.8cm), plantar/dorsi +/- 20°, toe-in/toe-out, swivel 30° left/right, and in height. Must choose size of Foot Holder in drop box.	X	X
Secure Foot Straps		Velcro® and D-ring adjustment holds feet securely in place at ankles & toes with multiple attachment slots for proper foot positioning. Straps are measured in length over top of foot, with approx. 1.5" (4cm) of Velcro® overlap. Must choose length of strap in drop box.	X	X
Multi-Adjustable Knee Pads		Multi-adjustable knee pads are adjustable in swivel and height +/75" (2cm), depth range of 5" (13cm) on size 1 models and 7" (18cm) on size 2 models. Must choose inside width of knee pad in drop box.	X	X

FEATURE	IMAGE	DESCRIPTION	MOD	MAX
Knee Strap	For Size 1 Frames:	1.5" (4cm) Velcro® strap with a D-Ring closure, adjustable in height. Size 1 fits circumference up to 15" (38cm). Size 2 fits circumference up to 19" (48cm) and has a 4"x6" (10x15cm) pad attached.	X	X
	For Size 2 Frames:			
Planar Pelvic Pad		Cushioned pad covered with Neoprene® provides pelvic support. Must choose pad size in drop box.	X	X

FEATURE	IMAGE	DESCRIPTION	MOD	MAX
Planar Upper Body Pad		Cushioned pad covered with Neoprene® provides upper body support. Must choose pad size in drop box	X	X
Hip Supports		Keeps the hips in symmetrical alignment. For Size 1: Pad size is 3.5"x3.5" (9x9cm), width is adjustable from 6"-10" (15-25cm). For Size 2: Pad size is 4.5"x4" (11x10cm), width is adjustable from 7"-13" (18-33cm) OR 12"-18" (30-46cm). Fore/aft adjustment of 4" (10cm). Must choose width adjustment range for Size 2 in drop box.	X	X
Padded Positioning Straps (2)		5"x8" (13x20cm) pad attached to a 2" (5cm) Velcro® belt with a D- Ring closure; fits circumference up to 24" (61cm) for Size1. 6"x10" (15x25cm) pad fitting circumference up to 35" (89cm) for Size 2.	X	X

FEATURE	IMAGE	DESCRIPTION	MOD	MAX
Lateral Supports	For Size 1 ONLY – Square Pads: For Size 1 &2 – Flat Pads:	Keeps the upper torso in symmetrical alignment. For Size 1 frame ONLY Square Pads are offered measuring 3.5"Hx3.5"W (9x9cm). For both Size 1 & 2 frames, Flat Pads measuring 4"Hx6.25"W (10x16cm) and Curved Pads measuring 4"Hx6"W (10x15cm) are offered. For Size 1 frame the style and range between pads (Narrow 4"-8" OR Wide 7"-11") must be chosen in drop down box. For Size 2 frames, choose style of pad in drop down box. Distance between laterals will be determined by your choice of Upper Body Support Planar Pad. Adjustable range of 7" - 12" (18 -30cm) with 7"W planar pads, 9" -14" (23 - 36cm) with 9"W planar pads.	X	X
	For Size 1 & 2 — Curved Pads:			

FEATURE	IMAGE	DESCRIPTION	MOD	MAX
Prone Face Aperture		Face Aperture has independently adjustable pads and supports the head. Pad size is 9"Hx3"W (23x8cm).	X	X
Multi-Adjustable Black Molded Tray		Tray provides anterior support and is adjustable in depth, height and angle. For Size 1 frames, Tray size is 11"Lx19"W (28x48cm). For Size 2 frames, Tray size is 21"Lx24"W (53x61cm).	X	X
Multi-Adjustable Clear Tray		Tray provides anterior support. For Size 1 frames: Tray size is 11"Lx19"W (28x48cm). For Size 2 Frames Tray size is 21"Lx24"W (53x61cm). Trays are adjustable in depth, height and angle.		X
Forearm Wings (On Size 1 Frame Only)		Forearm wings give additional area to rest arms when in standing position for size 1 frames only.	X	X

FEATURE	IMAGE	DESCRIPTION	MOD	MAX
Forearm Wings (On Size 1 Frame Only)		Forearm wings give additional area to rest arms when in standing position for size 1 frames only.		X
Elbow Pad (On Size 2 Frame Only)		Additional padded support surface on tray for size 2 frames only. Overall width 21.5" (55cm), pad extends back 3.75" (10cm) posterior from current tray edge for additional upper body/extremity support. Must choose size of cutout width in drop down box.	X	X
Push Handle for Tray	Size 1 Frame Size 2 Frame	Makes pushing the stander easier for the caregiver and aids in lifting and lowering the child.	X	X

FEATURE	IMAGE	DESCRIPTION	MOD	MAX
Foot Operated Gas Spring		The Gas Spring Lift Mechanism is used to adjust the angle and the position of the stander and is controlled by the Foot Pedal.	X	
Dual Control - Hand and Foot Operated Gas Spring		The Gas Spring Lift Mechanism is used to adjust the angle and the position of the stander. Dual Control enables position changes by foot pedal and hand control.		X
Gas Spring Lift Lockout		Locks gas spring in any position from lying down to fully upright and prone positions. Lift Lockout prevents any inadvertent position changes which can be useful if other children are present during standing sessions.		X

FEATURE	IMAGE	DESCRIPTION	MOD	MAX
Mast Receiver with Leg Abduction		Mast is the base for support and positioning components and allows for 0-30° of independently adjustable leg abduction.	X	X
Zing Swivel Locking Castors		Locking casters allow for easier steering and movement of the stander.	X	
Swivel Locking Castors with one Directional Locking Castor		Locking casters allow for easier steering and movement of the stander. The Directional Locking Castor locks the direction of one castor while leaving the other three free to swivel. This feature is very helpful when transporting the stander.		X